

Micro annular gear pumps with brushless motor

Programmable controller S-BL

For discrete dosage and continuous delivery



- High quality pump controller for continuous delivery and discrete dosage
- For mzl-pumps with special brushless motor
- Programming of controller with Windows® software »Motion Manager«
- Potentiometer for speed set
- Analog input 0-10 V
- 1 digital input, input is equipped with a switch
- 1 digital output, optionally programmable as input
- Two colored LED status indicator
- EEPROM memory
- RS-232 interface

The controller S-BL is recommended for high requirements in control of discrete and continuous delivery tasks in combination with a mzl-pump with special brushless motor.

The 16-bit microcontroller allows speed and position control for highly accurate dosage. The compact design on a PCB offers flexible installation. Process control link can be established via a

RS-232 interface. Motor speed or flow rate can be set either by an analog input (0-10 V) or a potentiometer mounted on the PCB. Programs for dosage can be saved in the memory.

Technical data

Control	PI-controller, speed and position control
Supply voltage U_B	24 V DC (12 – 30 V)
Speed	1 – 6000 rpm
Power	DIN 45323 socket, screw terminal
Pump connector	screw terminal, 8-pole
Serial interface	RS-232, SUB-D plug, 9-pole
Input # 1 (speed)	0 – 10 V
Error output (input # 2)	Open collector max. $U_B / 30 \text{ mA}$ no error: connected to GND as input: low 0...0.5 V / high 4 V... U_B
Digital inputs # 3	low 0...0.5 V / high 4...30 V input # 3 with switch
Program memory	6,600 Bytes
Protective class	IP 20
Dimensions (L x W x H)	approx. 112 x 85 x 36 mm
Weight	approx. 170 g

Subject to technical changes.

Contact

HNP Mikrosysteme GmbH
Bleicherufer 25 · D-19053 Schwerin

phone +49 385 52190-301
fax +49 385 52190-333

e-mail info@hnp-mikrosysteme.de
<http://www.hnp-mikrosysteme.de>

Software



All motor parameters for pump control can be set and saved with »Motion Manager« operating under Windows®. The program language is ASCII-based. Dosage programs can easily be typed on a computer and transferred to the EEPROM.

Several sample programs are supplied such as the triggering of different dosage programs with the internal PCB switches.

Configuration

```
SOR0      ;RS-232
SOR1      ;potentiometer
          ;or 0-10 V

V3000    ;speed 3000 rpm
V0        ;stop

LR10000   ;load 10 revolutions
M         ;start positioning

SP6000    ;set maximum speed
AC500     ;set acceleration
LPC400    ;load peak current
LCC200    ;load cont. current

GV        ;get command velocity
=> 3000   ;get actual velocity
=> 2998   ;get actual current
=> 200    ;get actual position
=> 1000   ;get configuration

EEPSAVE  ;save configuration
```

Set interface

Flow rate set via speed control

Dosage

Set motor parameters

Retrieve parameters and motor data


Save parameters to EEPROM

Programming

```
PROGSEQ  ;start program
A1        ;label 1
V1000    ;speed 1000 rpm
DELAY100  ;wait 1 sec
NP        ;notify position
LR 10000  ;load 10 revolutions
M         ;start motion
JMP1     ;jump to label 1
END      ;end program
```

Dosage program

Graphic online analysis



Sample chart: speed and motor current

Mode of operation

CONTMOD	standard mode
APCMOD	analog position control mode
STEPMOD	stepper motor mode

Selection of the control mode

Item number

66 02 01 05

controller S-BL for mzl-pumps with special brushless motor, null-modem cable, software »Motion Manager«, sample programs

Accessories

Power supply

external power supply 24 V DC, input voltage 100–240 V AC 50/60 Hz, with connector for controller S-BL

Multiplexer module

simultaneous operation of up to 255 pumps with a single RS-232 interface